
English

- PRODUCT INFORMATION:**
- 1) Barrel, Plunger & Tip Cap - Polypropylene
 - 2) Gasket - Bromo butyl rubber (Latex-free)
 - 3) Gasket Mixing – Silicon
 - 4) Barrel spray - Silicon (Medical grade)
 - 5) Sodium Chloride - 0.9% w/v sodium chloride Normal (saline) solution

INTENDED USE:

Easy Flush® Syringe is intended for flushing compatible intravenous administration sets and indwelling intravenous access devices.

CONTRAINDICATIONS:

- Do not use if the patient has a known allergy to any of its components, materials, or 0.9% w/v sodium chloride solution. It may lead to an allergic response resulting in anaphylaxis.
- Do not use if the intravenous administration sets and/or indwelling intravenous access devices are contaminated by microorganisms or the human body (skin) around the access device is inflamed
- Do not use in patients suffering from hypernatremia and fluid retention when the administration of sodium or chloride could be clinically detrimental.

PRECAUTIONS:

- EasyFlush® Syringe is designed to be used with ISO luer compliant components for intravenous applications. When handling the product care should be taken to avoid damage from handling.
- Check the product for visual defects/cracks.
- Gradation is provided on the syringe to inject the required amount of saline.
- Avoid crushing and crimping damage to the product due to the application of surgical instruments such as forceps.
- Store at a controlled temperature (10-40°C) and humidity (55% - 75%). Do not leave at freezing temperature.
- Check with drug manufacturer instructions for use to ensure compatibility with 0.9% w/v sodium chloride solution before use. (0.9% w/v sodium chloride solution is not compatible, follow the drug manufacturer's instructions for flushing practices, or first flush the vascular access device (VAD) with a compatible solution to remove traces of the medication in accordance with manufacturer and institution policies.
- Clinicians should consider the patient's specific medical conditions, treatment needs, age, and weight that may require restricted sodium or fluid intake when deciding to flush with 0.9% w/v sodium chloride injection.
- Saline flushes should be taken into account when prescribing fluids to not exceed fluid intake guidelines. Follow instructions for use, and institution procedures for flush administration

CAUTION:

- Carefully read all instructions prior to use.
- Do not use if the package is open or damaged.
- For single use only. Discard any partially used product.
- The solution and fluid path are sterile and non-pyrogenic if the tip cap is in place, the syringe is intact, and there is no evidence of leakage. Use proper aseptic technique

EXPIRATION PERIOD: The device has a shelf life of 3 years from the date of manufacturing. The product should not be used after the indicated sterility expiration date.

STORAGE:

- Store at a controlled temperature (10-40°C). Do not store at freezing temperature

INSTRUCTIONS FOR USE:

- Inspect the plastic wrapping. Do not use it if the plastic wrapping is damaged or not intact (Figure 1).
- Remove the plastic packaging by tearing along the perforation (Figure 2).
- Inspect the solution for any discoloration, cloudiness, haziness, or precipitate, or if the syringe is damaged (Figure 3).
- Ensure the cap is in place, then depress the plunger to release the stopper seal and activate the syringe. Improper activation may make the syringe difficult to use and/or cause plunger rod separation. Never draw back the plunger rod, as the product may become contaminated (Figure 4).
- Remove the tip cap by twisting it off using aseptic technique. Hold the syringe upright and prime it to expel any air bubbles if present (Figure 5).
- The syringe is now ready to use. Per institutional protocol, attach the flush syringe to the access device and flush (Figure 6).
- Inject the required amount of saline by pressing the syringe plunger slowly to avoid excessive pressure. If you experience plunger resistance, do not apply excessive force. Do not pull or bend the plunger sideways (Figure 7).

- Use according to the intravenous tubing or indwelling device manufacturer's recommendations. After flushing, gently remove the syringe from the access device (Figure 8).
- Discard the empty unit after use. Discard any unused portion. Do not reuse disposable syringes (Figure 9).

INTENDED PATIENT POPULATION:

EasyFlush® Syringe is intended to be used in patients irrespective of age and gender.

INTENDED USER:

EasyFlush® Syringe is intended to be used by trained healthcare professionals including doctors, registered practitioners, paramedics or nursing staff experienced in the use of vascular access devices.

HOW SUPPLIED:

EasyFlush® Syringe is supplied sterile and packed in a ribbon or paper pouch. Thirty (30) such pouches are enclosed within a printed duplex box. The overall packing contains eight (8) duplex boxes and external identification stickers.

WARNINGS:

- Do not reuse. Re-use may lead to infection or other illness.
- Do not re-sterilize the device by any method. Do not reprocess. It is not intended to be reprocessed and used again, even on the same patient.
- Reprocessing may compromise sterility, biocompatibility, and functional integrity of the device.
- SPM Medicare assumes no liability for the devices that have been re-use or re-sterilized.
- Do not use it if the product has been left at freezing temperature.
- Do not use the device after expiry date.
- Do not use it if the unit package or content is damaged.
- Store at specified temperature and humidity.
- Small parts are a potential choking hazard. After use, discard small parts according to the hospital protocol. Do not use it if the solution is cloudy or colored, contains a precipitate, or has any type of suspended particulate matter

EXPLANATION OF SYMBOLS:

	Caution		Medical Device		Non-Pyrogenic
	Manufacturer		Unique Device Identifier		Sterilized using irradiation
	Country of manufacturer		Single Sterile Barrier System with protective packaging outside		Handle with care
	Do not re-sterilize		Authorized representative in the European Community		Maximum stacking
	Consult instructions for use		Biological risks		Prescription Use
	Keep away from sunlight		Catalogue Number		Model number
	Do not re-use		Batch code		Single Sterile Barrier System
	Fluid Path		Use-by date		Date of Manufacture
	This side up		Temperature Limit		Fragile, handle with care
	Do not use hook		Do not use if package is damaged and consult instructions for use		Keep dry

Manufactured by:
SPM MEDICARE PVT. LTD.
(AN ISO 13485 Certified Company)
B-40, Phase-II, Noida, Gautam Buddha Nagar, Uttar Pradesh-201305, INDIA
Customer Care No. - +91-8407010118
E-mail: info@spmmedicare.com, Website: www.spmmedicare.com
Mfg. Lic. No. - MFG/MD/2023/000234
SIL by: IAS Shivram Institute of Industrial Research
19, University Road, Delhi-110007, India

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