



INSTRUCTIONS FOR USE

Oral Syringe

INSTRUCTIONS FOR USE

SPMM/RA/IND/42/ORLS/IFU, Rev-00, Date: 29/06/2026

Device Description

Oral Syringe with Tip Cap is used for the transport and administration of non-injectable liquid medication. The Tip Cap is used to cap the syringe to prevent leakage during transportation.



Intended Use

The Oral Dispensing Syringe are used for the transport and administration of non-injectable liquid medication.

Intended User

- Pharmacists, Nurses and Physicians experienced in the use of the syringe.
- Experience levels will be from novice to expert.

Instructions For Use

1. Accurately measure prescribed dosage using the top edge of the stopper. Syringe scales are indicated in milliliters (mL).
2. After filling syringe, attach tip cap which is provided separately in small bag.
3. Before dispensing, REMOVE TIP CAP and position syringe tip in the pocket between the cheek and the back teeth (molars).
4. To avoid choking, slowly dispense the medication by pushing gently on the plunger rod.
5. Discard syringe properly after use.

Warnings

INSTRUCTIONS FOR USE

SPMM/RA/IND/42/ORLS/IFU, Rev-00, Date: 29/06/2026

- Oral syringes are non-sterile product, and not intended to be sterilized prior to use.
- Do not use if the product or package has visible particulates.
- The Oral Syringe is to be used for Oral Use only.
- The Oral Syringe does not have any mating function except with Oral Tip Caps and is specifically incompatible with luer (6% taper) connectors and intravenous devices. When using tip cap, ensure secure connection to prevent leakage. Do not use if product is damaged.
- Potential Procedure Related Complications: Exposure to toxic medication.

Precautions

- Pharmaceutical products prepared into these syringes should be administered as soon as possible. For light sensitive medication, use amber colored oral syringe.
- Re-use may lead to infection or other illness/injury.

Disposal

Discard syringe properly after use according to applicable local waste disposal requirements.

LEGAL MANUFACTURER AND MANUFACTURING FACILITY



Manufacturer:

SPM Medicare Pvt, Ltd

B-40, Phase-II, Noida, Gautam Buddha Nagar, UP-201305, India.

Manufacturing Site 1- B-40, Phase-II, Noida,

Gautam Buddha Nagar, UP-201305, India.

Customer care: No.: +91-840-7070-718

Email: info@spmmedicare.com

Web Site: www.spmmedicare.com

EXPLANATION OF SYMBOLS-

Symbol	Explanation	Symbol	Explanation
	Caution		Catalogue Number
	Manufacturer		Batch code

INSTRUCTIONS FOR USE

SPMM/RA/IND/42/ORLS/IFU, Rev-00, Date: 29/06/2026

	Date & Country of manufacturer		Use-by date
	Medical Device		Do not use hook
	Consult instructions for use		Do not use if package is damaged
	Keep away from sunlight		Recycling
	Keep dry		Handle with care
	This side up		Model number

Revision History

Rev. No.	Date	Description of Change
00	29/06/2026	Initial release of Instructions for Use for Oral Syringe