

INSTRUCTIONS FOR USE
SPMM/RA/IND/52/SC/IFU, Rev. 00. Date: 29/06/2026

INSTRUCTION FOR USE

SUCTION CATHETER

 DATE: 29/06/2026	 DATE: 29/06/2026	 DATE: 29/06/2026
Prepared By Executive-	Reviewed By Head	Approved By Director/Auth. Signatory

INSTRUCTIONS FOR USE

SPMM/RA/IND/52/SC/IFU, Rev. 00. Date: 29/06/2026

Product Name:

Suction Catheter

Product Description:

The Suction Catheter is a sterile, single-use flexible tube designed for aspiration of secretions from the respiratory tract. The catheter consists of a soft atraumatic distal tip, suction tubing, color-coded connector, and control port. The device is inserted through the nose, mouth, endotracheal tube, or tracheostomy tube and connected to a vacuum source for suctioning procedures.

The catheter is designed to provide efficient secretion removal while minimizing trauma to the airway mucosa.

Intended Use:

The Suction Catheter is intended for the removal of respiratory secretions, mucus, blood, and other fluids from the airway and tracheobronchial tree of patients to maintain airway patency and facilitate respiration. The device is used in hospitals, intensive care units, emergency care, and home healthcare settings under professional supervision.

Shelf Life:

5 Years

Material of Construction:

- Medical Grade Polyvinyl Chloride (PVC)
- Polypropylene (PP)
- Polyethylene (PE)
- Medical Grade Colorant

Storage Condition

Store at controlled room temperature in a cool, dry, and clean place away from direct sunlight, excessive heat, moisture, and corrosive environments.

Principle of Operation

The catheter is connected to a vacuum source and inserted into the airway. Negative pressure generated by the suction apparatus removes secretions, mucus, blood, or fluids through the catheter lumen. The aspirated material is transported through the catheter and collected in an external suction collection system.

INSTRUCTIONS FOR USE

SPMM/RA/IND/52/SC/IFU, Rev. 00. Date: 29/06/2026

LEGAL MANUFACTURER AND MANUFACTURING FACILITY



Manufacturer:

SPM Medicare Pvt, Ltd

B-40, Phase-II, Noida, Gautam Buddha Nagar, Uttar Pradesh-201305, India.

Site 1- B-40, Phase-II, Noida, Gautam Buddha Nagar, Uttar Pradesh-201305, India.

Site 2- Industrial Building Plot no 26, Sector 158, Noida, Gautam buddha, Nagar, Uttar Pradesh, 201305.

Customer care: No.: +91-840-7070-718

Email: info@spmmedicare.com

Web Site: www.spmmedicare.com

EXPLANATION OF SYMBOLS

Symbol	Explanation	Symbol	Explanation
	Caution		Catalogue Number
	Manufacturer		Batch code
	Sterilized using Ethylene Oxide		Use-by date
	Do not re-sterilize		This side up
	Consult instructions for use		Handle with care
	Keep away from sunlight		Model number
	Do not re-use		Do not use hook
	Medical Device		Keep dry

INSTRUCTIONS FOR USE

SPMM/RA/IND/52/SC/IFU, Rev. 00. Date: 29/06/2026

Symbol	Explanation	Symbol	Explanation
	Model number		Do not use if package is damaged and consult <i>instructions for use</i>

Rev No	Description	Date
00	Initial Release of IFU	29/06/2026