



INSTRUCTIONS FOR USE

Sterile Latex Surgical Gloves (Powder Free)

INSTRUCTIONS FOR USE

SPMM/RA/IND/73/SG/IFU, Rev. No.: 00, Date: 11/09/2026

Product Name- Sterile Latex Surgical Gloves Powder Free

Classification- Class B

Intended Use

Sterile Latex Surgical Glove intended to protect the healthcare workers and patients from cross contamination when used in medical surgery.

Sterility – EO Sterilized/Gamma Sterilized

Caution

- This product contains natural rubber latex which may cause allergic reactions in some individuals. In case of latex allergy discontinue use and consult a physician.
- Sterile until package is opened or damaged.
- This glove does not protect the user from mechanical injuries, cold objects, heat or fire. It is also not intended for protection against electric shock, and radiation.

Warnings

- The performance levels obtained in permeation testing does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves
- Before usage, inspect the gloves for any defect or imperfections
- For single use only

INSTRUCTIONS FOR USE

SPMM/RA/IND/73/SG/IFU, Rev. No.: 00, Date: 11/09/2026

Directions For Use

1. Dry hands, select size, and check for damage.
2. Open package and put on gloves aseptically. Replace if damaged or contaminated.
3. Single use only. Do not reuse or re-sterilize. Dispose after use.

Storage instruction

Storage in cool, dry place away from direct sunlight, in the temperature of range 5-30⁰ C
Glove characteristics will not be changed significantly on storage, when it is stored as per recommended storage conditions.

EXPLANATION OF SYMBOLS:

Symbol	Explanation	Symbol	Explanation
	Caution		Catalogue Number
	Manufacturer		Batch code
	Date & Country of manufacturer		Use-by date
	Do not re-use		Medical Device
	Consult instructions for use		Do not use if package is damaged
	Keep away from sunlight		Keep dry
	Do not re-sterilize		

Revision History:

Rev. No	Reason for Change	Date
00	Initial Release of IFU	11/09/2026