

Committed to Care



Surgical Disposable

DISPOSABLE TROCAR NEEDLE


Easy-Trocar Suprapubic Needle Surgical

DISPOSABLE OPTICAL TROCAR NEEDLE

BLADELESS



TECHNICAL SPECIFICATIONS



Features :

- Suprapubic needle surgical used in minimally invasive surgery.
- Optical entry for endoscope.
- Smooth entry with low insertion force into the abdominal cavity.
- Funnel design enables smooth instrument exchange and specimen removal.
- Clear housing provides surgeon precise visibility of shield position.
- Thread design maintains the seal & helps ensure a leak free operation.
- The ergonomic design of the trocar thread not only makes instruments stable but also reduces harm to skin.
- Pattern and color design is pleasing and comfortable to eyes.

Technical Specification :

Brand	Size (mm)	Product code	Specification	Qty. Inner /Outer
EASY TROCAR Suprapubic Needle Surgical (DISPOSABLE OPTICAL TROCAR NEEDLE)	5x100	101ABSP1	Suprapubic Needle Surgical, paper pouch pack, size : 5x100mm	10/40
	10x100	101ABSP2	Suprapubic Needle, paper pouch pack, size : 10x100mm	10/40
	12x100	101ABSP3	Suprapubic Needle, paper pouch pack, size : 12x100mm	10/40
	5x100	101ABSB1	Suprapubic Needle, Individual inner box pack : 5x100mm	1/45
	10x100	101ABSB2	Suprapubic Needle, Individual inner box pack : 10x100mm	1/45
	12x100	101ABSB3	Suprapubic Needle, Individual inner box pack : 12x100mm	1/45

SUPRAPUBIC NEEDLE SURGICAL (Disposable Trocar Needle)



Individual Pack Inner box

Easy-Cision

Soft Tissue Trephine

DISPOSABLE CIRCUMCISION STAPLER

DOUBLE HANDLE



Application :

- Used for male circumcision in urology


Features :

- Disposable sterilizing product. Prevent various infections, no complications after operation.
- Medical staplers are designed to be user-friendly for healthcare professionals. & provide precise and consistent closure of tissues to promote proper healing.
- The incision margins are neat, and the postoperative recovery is fast.
- The operation is standardized, the instruments are simple to operate and the procedure is fast
- Accurate control of the scope of excision, regular cutting edge, beautiful, small scar, fast recovery.
- Different types of rings are available according to individual differences.
- Individual blister pack & Outer Pack : 1Pcs x 20Boxes =20 Pcs

C-type Foreskin Stapler Parameter Table

S.No.	Product Code	Size / Specification	Cutting knife Diameter	Tol. ±	Nail height	Tol. ±	Suture Needle Button
1	102DB12B	Ø 12	13	±2	3.0	±5	10
2	102DB15B	Ø 15	15				11
3	102DB18B	Ø 18	18				13
4	102DB22B	Ø 22	23				16
5	102DB26B	Ø 26	25				18
6	102DB30B	Ø 30	28.5				20
7	102DB32B	Ø 32	32				22
8	102DB34B	Ø 34	34				24

Easy-Clip Titanium Ligating Clips



Application :

- A surgical device designed to apply a hemostatic clip for ligation to aid in performing.

Features :

- Safe locking device with a more secure closure
- No supporting bacterial growth
- Color-coded clips and corresponding color-coded applicators for different sizes of clips
- For use with open and endoscopic single clip applicators.
- Features proprietary lateral and transverse grooves designed for secure fixation on the structure

Technical Specification :

Brand	Product Code	Color	Specification	Package
Easy-Clip Titanium Ligating Clips	103EAB	Yellow	Titanium Ligating Clip Size -100 (Small)	6 Clips/Cartridge 18 Cartridge/Box
	103EAB	Blue	Titanium Ligating Clip Size -200 (Medium)	6 Clips/Cartridge 18 Cartridge/Box
	103EAB	Green	Titanium Ligating Clip Size -300 (Medium Large)	6 Clips/Cartridge 18 Cartridge/Box
	103EAB	Orange	Titanium Ligating Clip Size -400 (Large)	6 Clips/Cartridge 18 Cartridge/Box

Our Certificates

ISO Certificate

ITC
Management Systems Certification Body
Institute for Testing and Certification, Inc.
Vřida Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic
www.itc299.cz

CERTIFICATE
No. 21 0097 SJ

We confirm on the basis of a performed audit that company
SPM Medicare Pvt. Ltd.
B-40, Phase II, Noida, Gautam Buddha Nagar, Uttar Pradesh - 201305, India
Company Reg. No.: U33110DL2014PTC269667

has implemented and documented a functional quality management system
in compliance with the requirements of the standard

EN ISO 13485:2016

Covering the following activities:
Manufacturing, Sales and Distribution of Single Use Sterile Medical Devices for
Injection, Infusion, Transfusion & Medical Apparatus

The Certificate is issued on the basis of the results mentioned in Audit Report No.
23340034/2021. The Certificate validity is conditioned by positive results of surveillance audits,
which the certified company committed to undergo.
During use of the Certificate the Certificate Holder undertakes to follow the Rules of Use of the Certificate.
This document is publicly available on www.itc299.cz

Ing. Pavel Vaněk
Head of Certification Body

Date of Issue: 20.12.2021
Valid until: 19.12.2024

License Certificate

FORM MD-5
[See sub-rule (4) of rule 20 and sub-rule (6) of rule 20]

Licence to Manufacture for Sale or for Distribution of Class A or Class B medical device

Licence Number: MFG/MD/2020/000315
Endorsement No. 10

1. M/s SPM Medicare Pvt. Ltd. B-40, Phase 2, Noida, Gautam Buddha Nagar, Uttar Pradesh (India) - 201305
Telephone No.: 9711008401 FAX: 120 4634949 has been licensed to manufacture for sale or for distribution the below
listed medical device(s) at the premises situated at M/s SPM MEDICARE, B-40, PHASE-2nd, Noida, Uttar Pradesh,
Gautam Buddha Nagar, Uttar Pradesh (India) - 201305 Telephone No.: 9711008401 FAX: 01204270212

2. Details of medical device(s) [Annexed]

3. This licence is subject to the provisions of the Medical Devices Rules, 2017 and conditions prescribed therein.

ANNEXURE

S.No.	Details Of Device(s)
1	Generic Name: Soft Tissue Trephine (Disposable Circumcision Stapler) Model No.: Soft Tissue Trephine (Disposable Circumcision Stapler) - Single Handle Support - 12mm, 15mm, 18mm, 22mm, 26mm, 30mm, 32mm, 34mm, 36mm, Soft Tissue Trephine (Disposable Circumcision Stapler) - Double Handle Support - 12mm, 15mm, 18mm, 22mm, 26mm, 30mm, 32mm, 34mm, 36mm Intended Use: A cylindrical or conical saw used to resect disc of tissues other than bones. It comes in various sizes and configurations according to the size and hardness of the tissue to be resected. Class of medical device: Class B Material of construction: Stainless Steel, Silicon, 300-GSM Paper, Non-Woven Dimension (if any): 12-36 mm Shelf life: 5 Year Sterile or Non-sterile Sterilized Brand Name (if registered under the Trade Marks Act, 1999): NA

Place: Lucknow
Date: 08-Dec-23

Shashi Mohan
Gupta
(Shashi Mohan Gupta)
Drug Licensing and Controlling Authority
Uttar Pradesh
State Licensing Authority

Digitally signed by Shashi Mohan
Gupta
Date: 2023.12.08 16:38:43 +05'30'

License Certificate

FORM MD-5
[See sub-rule (4) of rule 20 and sub-rule (6) of rule 20]

Licence to Manufacture for Sale or for Distribution of Class A or Class B medical device

Licence Number: MFG/MD/2020/000315
Endorsement No. 10

1. M/s SPM Medicare Pvt. Ltd. B-40, Phase 2, Noida, Gautam Buddha Nagar, Uttar Pradesh (India) - 201305
Telephone No.: 9711008401 FAX: 120 4634949 has been licensed to manufacture for sale or for distribution the below
listed medical device(s) at the premises situated at M/s SPM MEDICARE, B-40, PHASE-2nd, Noida, Uttar Pradesh,
Gautam Buddha Nagar, Uttar Pradesh (India) - 201305 Telephone No.: 9711008401 FAX: 01204270212

2. Details of medical device(s) [Annexed]

3. This licence is subject to the provisions of the Medical Devices Rules, 2017 and conditions prescribed therein.

ANNEXURE

S.No.	Details Of Device(s)
1	Generic Name: Soft Tissue Trephine (Disposable Circumcision Stapler) Model No.: Soft Tissue Trephine (Disposable Circumcision Stapler) - Single Handle Support - 12mm, 15mm, 18mm, 22mm, 26mm, 30mm, 32mm, 34mm, 36mm, Soft Tissue Trephine (Disposable Circumcision Stapler) - Double Handle Support - 12mm, 15mm, 18mm, 22mm, 26mm, 30mm, 32mm, 34mm, 36mm Intended Use: A cylindrical or conical saw used to resect disc of tissues other than bones. It comes in various sizes and configurations according to the size and hardness of the tissue to be resected. Class of medical device: Class B Material of construction: Stainless Steel, Silicon, 300-GSM Paper, Non-Woven Dimension (if any): 12-36 mm Shelf life: 5 Year Sterile or Non-sterile Sterilized Brand Name (if registered under the Trade Marks Act, 1999): NA

Place: Lucknow
Date: 08-Dec-23

Shashi Mohan
Gupta
(Shashi Mohan Gupta)
Drug Licensing and Controlling Authority
Uttar Pradesh
State Licensing Authority

Digitally signed by Shashi Mohan
Gupta
Date: 2023.12.08 16:38:43 +05'30'



SPM Medicare Pvt. Ltd.

(AN EN ISO 13485 Certified Company)

Corporate & Manufacturing Address:

B-40, Phase-II, NOIDA, Gautam Buddha Nagar - 201305 U.P., INDIA

Customer Care: +91-8407070718 E-mail: info@spmmedicare.com

Website: www.spmmedicare.com