


Committed to Care



Surgical Disposable


TROCAR (DISPOSABLE)

Easy-Trocar Trocar (Disposable)

BLADELESS



TECHNICAL SPECIFICATIONS



Features :

- Trocar surgical used in minimally invasive surgery.
- Optical entry for endoscope.
- Smooth entry with low insertion force into the abdominal cavity.
- Funnel design enables smooth instrument exchange and specimen removal.
- Clear housing provides surgeon precise visibility of shield position.
- Thread design maintains the seal & helps ensure a leak free operation.
- The ergonomic design of the trocar thread not only makes instruments stable but also reduces harm to skin.
- Pattern and color design is pleasing and comfortable to eyes.

Technical Specification :

Brand	Size (mm)	Product code	Specification	Qty. Inner /Outer
EASY TROCAR	5x100	101ABSP1	Paper pouch pack, size : 5x100mm	10/40
	10x100	101ABSP2	Paper pouch pack, size : 10x100mm	10/40
	12x100	101ABSP3	Paper pouch pack, size : 12x100mm	10/40
TROCAR (DISPOSABLE)	5x100	101ABSB1	Hard Blister, Individual inner box pack Size : 5x100mm	1/45
	10x100	101ABSB2	Hard Blister, Individual inner box pack Size : 10x100mm	1/45
	12x100	101ABSB3	Hard Blister, Individual inner box pack Size : 12x100mm	1/45

TROCAR (DISPOSABLE)



Individual Pack Inner box

CIRCUMCISION (DISPOSABLE)

Easy-Cision

Soft Tissue Trephine

DISPOSABLE CIRCUMCISION STAPLER



DOUBLE HANDLE



Single row with Multiple Pin



Diffirent Cup size

Application :

- Used for male circumcision in urology

Features :

- Disposable sterilizing product. Prevent various infections, no complications after operation.
- Medical staplers are designed to be user-friendly for healthcare professionals. & provide precise and consistent closure of tissues to promote proper healing.
- The incision margins are neat, and the postoperative recovery is fast.
- The operation is standardized, the instruments are simple to operate and the procedure is fast
- Accurate control of the scope of excision, regular cutting edge, beautiful, small scar, fast recovery.
- Different types of rings are available according to individual differences.
- Individual blister pack & Outer Pack : 1Pcs x 20Boxes =20 Pcs



C-type Foreskin Stapler Parameter Table

S.No.	Product Code	Size / Specification	Cutting knife Diameter	Tol. ±	Nail height	Tol. ±	Suture Needle Button
1	102DB12B	Ø 12	13	± 2	3.0	± 5	10
2	102DB15B	Ø 15	15				11
3	102DB18B	Ø 18	18				13
4	102DB22B	Ø 22	23				16
5	102DB26B	Ø 26	25				18
6	102DB30B	Ø 30	28.5				20
7	102DB32B	Ø 32	32				22
8	102DB34B	Ø 34	34				24



SPM Medicare Pvt. Ltd.

(AN EN ISO 13485 Certified Company)

Corporate & Manufacturing Address:

B-40, Phase-II, NOIDA, Gautam Buddha Nagar - 201305 U.P., INDIA

Customer Care: +91-8407070718 E-mail : info@spmmedicare.com

Website : www.spmmedicare.com